

# CTA Chest/abd/pelvis 64 Toshiba

<b>Indications</b>	trauma, acute aortic syndrome, suspected aneurysm/dissection										
<b>Diagnostic Task</b>	Detect aneurysms, aortic dissections										
<b>Scan mode</b>	Helical										
<b>Position/Landmark</b>	Head or feet first-Supine 1cm superior to shoulder										
<b>Topogram</b>	AP mA50 kV120 /Lat mA 70 kV120										
<b>kVp/Reference mass</b>	135kv Sure Exp 3D(80-550)										
<b>Rotation time/pitch</b>	0.5\0.828										
<b>Detector Configuration</b>	64x0.5										
<b>Table Speed/Increment</b>	26.5										
<b>Dose reduction</b>	Sure Exp 3D										
<b>Allowed CTDI ranges*</b>	7mGy-50mGy										
<b>XR29 Dose Notification value</b>	50mGy										
<b>Helical Set #1</b>	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;"></td> <td style="width: 15%; text-align: center;">body</td> <td style="width: 15%; text-align: center;">thickness</td> <td style="width: 15%;"></td> <td style="width: 15%; text-align: center;">recon</td> </tr> <tr> <td style="text-align: center;">recon</td> <td style="text-align: center;">part</td> <td style="text-align: center;">spacing</td> <td style="text-align: center;">algorithm</td> <td style="text-align: center;">destination</td> </tr> </table>		body	thickness		recon	recon	part	spacing	algorithm	destination
		body	thickness		recon						
recon	part	spacing	algorithm	destination							
1	chest	1.5mmx 1.5mm	standard	pacs							
<b>Helical Set #2</b>	1	chest	2mmx 2mm	standard	pacs						
	2	lung	1mmx1mm	lung	pacs						
	3	sag cap	2mmx2mm	standard	pacs						
	4	coronal cap	2mmx2mm	standard	pacs						
	5	axial mip lung	10mmx2mm	standard	pacs						
	6	MIP coronal aorta	5mmx2mm	standard	pacs						
	6	MIP sag aorta	5mmx2mm	standard	pacs						
	8	thin cap	1mmx0.8mm	standard	pacs/TR						
<b>When super D or stereo chest</b>	** When super D or stereo chest ordered and on all FHS pulmonary patients**										
<b>Scan Start/end location</b>	NC 2cm superior to lung apices// arterial 2cm superior to lung apices NC through hepatic dome// lesser trochanters										
<b>DFOV</b>	40cm										
<b>IV contrast volume/type</b>	100ml isovue 370 3-4cc/sec										
	Performed as directed by the supervising radiologist										
<b>Scan delay</b>	Surestart										
	bolus tracking in the descending aorta(level just inferior to carina)										
	Comments: Being able to locate the descending aorta is important. The monitoring phase will not trigger properly and the scan will not start correctly if the roi is not placed on the correct anatomy.										
	<b>Approximate Values for CTDIvol</b>										
	Patient size	weight(kg)	weight(lbs)	CTDIvol(mGy)							
	SMALL	50-70	110-155	4-10							
	AVERAGE	70-90	155-200	8-16							
	LARGE	90-120	200-265	14-22							
<b>NOTE*</b>	*The AAPM recommended NEMA XR29 Dose Notification Value for an adult torso is 50mGy. Dose Notification levels less than the AAPM recommended can be set. The maximum CTDI vol should match the dose notification value. Exams with CTDI vol values less than the minimum										

allowed range should not be performed unless approved by a radiologist.

