

# ROUTINE NECK/CHEST wo 16 Emotion

Indications	For abdomen pain, lymphoma, restage ca, weight loss, fatigue,						
Diagnostic Task	Detect masses, free fluid, abscess, mets						
Scan mode	Helical						
Position/Landmark	2cm superior to xiphoid/Inspiration						
Topogram	AP 25mA 130kV Lat 25mA 130kV						
kVp/Reference mass	130kv 120mas-110kv if pt over 140lbs						
Rotation time/pitch	neck 0.8/0.75 C/A/P 0.6/0.8						
Detector Configuration	neck 16x1.2 C/A/P 16x1.2						
Table Speed/Increment	neck 14.4 C/A/P 15.36						
Dose reduction	CareDose 4D						
Allowed CTDI ranges*	7mGy-50mGy						
XR29 Dose Notification value	50mGy						
Helical Set#1	body		thickness		recon		
Chest	recon	part	spacing	kernel	window	destination	
arms up	1	chest	2mmx2mm	31medium	smooth	Mediastinum	pac
	2	lung	1.5mmx1.5mm	70sharp		lung	pac
	3	sag chest	2mmx2mm	31medium	smooth	mediastinum	pac
	4	coronal chest	2mmx2mm	31medium	smooth	mediastinum	pac
	7	axial mip lung	10mmx2mm	b20f	smooth	lung	pac
Helical Set#2	body		thickness		recon		
arms down	recon	part	spacing	kernel	window	destination	
	1	neck	2mmx 2mm	31medium	smooth	mediastinum	pac
	2	coronal neck	2mmx2mm	31medium	smooth	mediastinum	pac
	3	sag neck	2mmx2mm	31medium	smooth	mediastinum	pac
Scan start	Chest-1cm superior to shoulder/			neck-top of orbital roof			
End location	L1			/ neck base			
FOV	40cm			20cm			
	decrease appropriately						
IV contrast-split bolus	na						
NOTE*	na						
	MARK AREA OF PAIN WITH BB						
	Approximate Values for CTDIvol						
	Patient size	weight(kg)	weight(lbs)	CTDIvol(mGy)			
	SMALL	50-70	110-155	10-17			
	AVERAGE	70-90	155-200	15-25			
	LARGE	90-120	200-265	22-35			
NOTE*	*The AAPM recommended NEMA XR29 Dose Notification Value for an adult torso is 50mGy. Dose Notification levels less than the AAPM recommended can be set. The maximum CTDI vol should match the dose notification value. Exams with CTDI vol values less than the minimum allowed range should not be performed unless approved by a radiologist.						

