

NON CON ABDOMEN/PELVIS Siemens GO ALL

| | | | | | | |
|-------------------------------------|---|------------------|-------------|--------------|---------|-------------|
| Indications | For abdomen pain, vomiting, bloating | | | | | |
| Diagnostic Task | Detect diverticulitis, free fluid, appendicitis, obstruction | | | | | |
| Scan mode | Helical | | | | | |
| Position/Landmark | 2cm superior to xiphoid/Inspiration | | | | | |
| Topogram | AP 110kV 15mA | | | | | |
| kVp/Reference mass | 120kv 99mas | | | | | |
| Rotation time/pitch | 0.5/0.8 | | | | | |
| Detector Configuration | 32x0.7 | | | | | |
| Table Speed/Increment | 17.92 | | | | | |
| Dose reduction | CareDose 4D, SAFIRE 3 | | | | | |
| Allowed CTDI ranges* | 7mGy-50mGy | | | | | |
| XR29 Dose Notification value | 50mGy | | | | | |
| Helical Set #1 | recon | part | spacing | kernel | window | destination |
| | | body | thickness | | | recon |
| | recon | part | spacing | kernel | window | destination |
| | 1 | abd/pelvis axial | 2mmx 2mm | Br40 | abdomen | pacs |
| | 2 | abd/pelvis Cor | 2mmx 2mm | Br40 | abdomen | pacs |
| 3 | abd/pelvis Sag | 2mmx 2mm | Br40 | abdomen | pacs | |
| Scan Start/end location | 1cm superior to diaphragm | | | | | |
| | lesser trochanters | | | | | |
| DFOV | 40cm | | | | | |
| | decrease appropriately | | | | | |
| IV contrast volume/type | none | | | | | |
| Scan delay | scanned during valsalva if looking for hernia | | | | | |
| | WITH ORAL CONTRAST ONLY | | | | | |
| | Approximate Values for CTDIvol | | | | | |
| | Patient size | weight(kg) | weight(lbs) | CTDIvol(mGy) | | |
| | SMALL | 50-70 | 110-155 | 10-17 | | |
| | AVERAGE | 70-90 | 155-200 | 15-25 | | |
| | LARGE | 90-120 | 200-265 | 22-35 | | |
| NOTE* | <p>*The AAPM recommended NEMA XR29 Dose Notification Value for an adult torso is 50mGy. Dose Notification levels less than the AAPM recommended can be set. The maximum CTDI vol should match the dose notification value. Exams with CTDI vol values less than the minimum allowed range should not be performed unless approved by a radiologist.</p> | | | | | |

