

SINUS COMPLETE HELICAL SUMMARY

Indications	Trauma, Pain, Swelling, sinus pressure, sinus drainage, cough				
Diagnostic Task	Detect fractures, edema, masses, or sinus infections				
Scan mode	Helical				
Helical Set	body		thickness		recon
	recon	part	spacing	algorithm	destination
	1	Sinus	.625mmx .625mm	bone	pacs
	2	Sinus large FOV	1mmx 1mm	standard	pacs
	3	coronal sinus	1mmx1mm	bone	pacs
	4	sag sinus	1mmx1mm	bone	pacs
Scan Start/End	1cm inferior to maxillary sinus				
	1cm superior to frontal sinus				
DFOV	Bone 18cm/Standard <25cm				
Angle	none				
IV contrast volume/rate	80ml isovue 370, 2cc/sec if needed				
Scan Delay	60 seconds/Performed as directed by the supervising radiologist				
	Mark Right side of patient with BB.				