

Peds Scoliosis Spine

Reviewed By: Justin Waters, MD

Last Reviewed: January 2018

Contact: (866) 761-4200, Option 1

General: If post operative lumbar spine, do with contrast.

If for myelopathy or multiple sclerosis, do with contrast

Sagittal imaging must extend through the facet joints into distal transverse processes.

All axial imaging should be contained within one sequence. If the images need to be angled, please angle within a single axial set of images as opposed to generating multiple axial sets with different angles.

Sequences:

- Sag T2 of entire spine for counting purposes
- Sag T1 cervical
- Sag T1 thoracic
- Sag T1 lumbar
- Sag T2 cervical
- Sag T2 thoracic
- Sag T2 lumbar
- Cor T1 cervical
- Cor T1 thoracic
- Cor T1 lumbar
- Ax T1 through major curve
- Ax T2 through major curve