

TRA-MINW

Routine Pediatric Brain <1y w/Contrast

Reviewed By: Douglas Hartz, MD

Last Reviewed: January 2018

Contact: (866) 761-4200, Option 1

General Comments: Chiari I malformation - Consult with a radiologist about additional imaging:

Add a Sag CSF flow study at foramen magnum.

Add Sag T2 of entire spine for evaluation of syrinx.

Sequences:

- Sag MPR
- Ax Diffusion
- Ax FLAIR
- Ax T2*
- Ax T1
- Ax T1 IR
- Cor T2
- Ax T2 post
- Ax T1 post
- Sag MPR post
- Coronal recon MPR post
- Axial recon MPR post