



Multiple Sclerosis WO

Reviewed By: Justin Waters, MD

Last Reviewed: January 2020

Contact: (866) 761-4200, Option 1

Use: Confirmed multiple sclerosis only. If suspected then perform Routine Adult Brain MRI with addition of Sag FLAIR.

Sequences:

All sequences other than MPR do at 3 mm

- Ax Diffusion
- Ax T2
- Ax FLAIR
- Ax T1
- Sag FLAIR
- Sag MPR
- Coronal recon MPR
- Axial recon MPR