



Peds Screening Spine for Drop Mets (New Tumor)

Reviewed By: Justin Waters, MD

Last Reviewed: January 2020

Contact: (866) 761-4200, Option 1

General: If for myelopathy or multiple sclerosis, do with contrast

Sagittal imaging must extend through the facet joints into distal transverse processes.

All axial imaging should be contained within one sequence. If the images need to be angled, please angle within a single axial set of images as opposed to generating multiple axial sets with different angles.

Sequences:

- Sag T2 of entire spine for counting purposes
- Sag T1 post cervical
- Sag T1 post thoracic
- Sag T1 post lumbar
- Ax T1 post cervical
- Ax T1 post thoracic
- Ax T1 post lumbar