



Pituitary Mass w/Contrast

Reviewed By: Justin Waters, MD

Last Reviewed: January 2020

Contact: (866) 761-4200, Option 1

Use: Pituitary mass other than microadenoma

Sequences:

- Sag T1 (Sag MPR for peds)
- Ax Diffusion
- Ax FLAIR
- Sag T2 thin thru sella
- Cor T1 thin thru sella
- Sag T1 thin thru sella
- Ax T2 post
- Cor T1 post thin thru sella
- Sag T1 post thin thru sella
- Sag T1 MPR post thru brain
- Coronal recon MPR post
- Axial recon MPR post