



Trigeminal Neuralgia WO

Reviewed By: Justin Waters, MD

Last Reviewed: January 2020

Contact: (866) 761-4200, Option 1

FOV: Similar routine for IACs, but centered at Meckel's Cave

Only to be used if there is a strong contraindication for contrast.

Sequences:

- Sag T1
- Ax Diffusion
- Ax FLAIR
- Ax T2
- Ax CISS Meckel's cave
- Ax T1 thin Meckel's cave
- Cor T1 fat sat thru Meckel's cave and brainstem