

Mediastinal mass MR Chest Without and With

Reviewed By: Spencer Lake, MD
Last Reviewed: October 2022

Contact: (866) 761-4200, Option 1

NOTE: Call radiologist before starting patient if any questions re: protocol

Indications: Mediastinal mass

Oral contrast: None.

Coil: Phased array body coil or cardiac coil centered over the mid chest.

Coverage: Mediastinum

Intravenous contrast: Single dose gadolinium @ 2 cc / sec (MultiHance)

1. **3-plane localizers with breath hold**
2. **Axial T2W SSFSE**
 - a. ECG-gated
3. **Coronal T2W SSFSE**
 - a. ECG-gated
4. **Axial T2W STIR**
 - a. Multi-breath hold
 - b. ECG-gating preferred but not required
5. **Axial balanced SSFP**
 - a. ECG-gated
6. **Axial T1W dual-echo in- and opposed-phase**
 - a. 3 mm slice thickness
7. **Diffusion-weighted imaging**
 - a. ECG-gated
 - b. b=0, 100, and 800

8. Precontrast T1W fast GRE – axial only
9. Postcontrast T1W fast GRE
 - a. Axial acquisitions at 20-30 sec, 60-70 sec, 3 min, and 5 min (with subtraction sequences)
 - b. Sagittal and coronal acquisitions between 3 min and 5 min

Check images with radiologist if any concerns.
