



Cervical Spine WO

Reviewed By: Justin Waters, MD

Last Reviewed: January 2024

Contact: (866) 761-4200, Option 1

General: If for myelopathy or multiple sclerosis, do with contrast

Sagittal imaging must extend through the facet joints into distal transverse processes.

All axial imaging should be contained within one sequence. If the images need to be angled, please angle within a single axial set of images as opposed to generating multiple axial sets with different angles.

Coverage: Sagittal images from skull base through T1-T2 only

- Sag T2
- Sag T1
- Sag STIR
- Ax T2
- Ax T2 SPACE if scanner allows. If no, do Ax MEDIC. If neither is possible, do Ax T2 GRE

Clinical Indications:

*Ordered by specialist or non-specialist as MR C-Spine WO

- Degenerative disease
- Disc herniation
- Neck pain
- Arm/shoulder pain and/or weakness