



Internal Auditory Canals (IAC) w/Contrast

Reviewed By: Justin Waters, MD

Last Reviewed: October 2023

Contact: (866) 761-4200, Option 1

Use: Asymmetric hearing loss, R/O acoustic neuroma, routine Kevin Kennedy case

DO NOT use for dizziness, use routine brain.

Sequences:

- Sag T1
- Ax Diffusion
- Ax Non-EPI Diffusion posterior fossa (if available)
- Ax FLAIR
- Ax CISS posterior fossa
- Ax T1 thin posterior fossa
- Cor T1 fat sat thru IAC
- Ax T2 post
- Ax T1 post posterior fossa
- Cor T1 fat sat post thru IAC
- Sag MPR post thru brain

Clinical Indications:

*Ordered by specialist or non-specialist as MR Brain W & WO, attention IAC

- Hearing loss