



## Multiple Sclerosis WO

**Reviewed By:** Justin Waters, MD

**Last Reviewed:** January 2024

**Contact:** (866) 761-4200, Option 1

**Use:** Confirmed multiple sclerosis only. If suspected then perform Routine Adult Brain MRI with addition of Sag FLAIR fatsat.

### **Sequences:**

All sequences other than MPR do at 3 mm

- Ax Diffusion
- Ax T2
- Ax FLAIR
- Ax T1
- Sag FLAIR fatsat
- Sag MPR
- Coronal recon MPR
- Axial recon MPR